

CUSTOM POSITIONING BLOCKS




Advanced Positioners™ Made To Order For You

Advanced Positioners™ are the highest quality available in today's medical market.

AADCO Medical Inc. offers a wide variety of durable polyfoam positioning blocks. If you cannot find a positioner to suit your specific need, we offer the service of making custom designed positioners.

To Order:

Many options are available for foam filled material, covering material, seam types, etc. The list below represents most of the available options.



Seam Types:

1. Standard Sewn
2. Waterproof (heat sealed blue vinyl only)

Filler Materials:

1. High density polyfoam
2. Durafoam™ (must be laminated in thickness over 2½")
3. PolyFoam F/R™ foam, special fire retardant foam, ester type 2 lb. density

Covering Materials:

1. Blue nylon-waterproof-highly durable
2. Blue vinyl (waterproof seams)
3. Special bonded urethane coating (Durafoam™ only)
4. Beige vinyl
5. Conductive Covering
6. Coated PolyFoam (water resistant and cleanable)

ORDER NO. PB-999

Note: Custom products are quoted upon request and are not returnable. We encourage you to call to discuss your specific need and make full use of our technical assistance before ordering. Thank You

Order Form

Dimensions					
	A	B	C	D°	E°
1					
2					
3					

Please fill out form completely. If space is inadequate, attach an additional sheet of paper.

1. Application _____
2. Specific use (if any) _____
3. Filler _____
4. Covering (if any) _____
5. Seams _____
6. Department _____
7. Your Name _____
8. Telephone No. _____
9. Purchase Order No. (if used) _____

Customer signature