


# CUSTOM TABLE PADS



**O.R., X-ray, Angio, Cardiac Cath,  
Nuclear, Cat Scan, E.R., Ultra-Sound,  
MRI, Electrophysiology,  
"We Make Them All"**

Advanced Table Pads & Systems are the highest quality available in today's medical market. Our Table Pad Systems are a unique combination of straps & fasteners with custom table pads.

The AADCO Medical Surgical Catalog offers a wide variety of standard size table pads. If you cannot find a table pad to suit your specific need, we offer the service of making custom designed table pads and table pad systems.

Many options are available for seam types, filler material, covering material, safety straps and Velcro® type fasteners

## Order Form

Table Pad Dimensions (Rectangles only, otherwise diagram below)

Side	A	B	C
Measurement			

Please fill out form completely. If space is inadequate, draw a small diagram below or append an additional sheet of paper.

1. Facility/Department \_\_\_\_\_
2. Your Name \_\_\_\_\_
3. Telephone No. \_\_\_\_\_
4. Purchase Order No. (if used) \_\_\_\_\_
5. Table Type & Model No \_\_\_\_\_

6. Diagram

**Please use one copy of this form for each pad being ordered**

### 7. Material Choice

**Filler Materials:** (check one)

- Contour Foam™ Available in 1.5" and 2" thicknesses.
- Super-support™ foam (must be laminated in thicknesses over 1" up to 2")
- Pressure Displacement Foam™ Available in thicknesses from 1" to 8"
- Lami-Foam™ available in 2" and 2.5" thicknesses
- Soft-Support Foam™ available in 1.5", 2" and 2.5" thicknesses
- Lo-Sorb Foam™ available in 1", 1.5", 2" thicknesses

### 8. Upholstery Choice

**Covering Materials:** (check one)

- Black nylon-waterproof - highly durable (water resistant seams)
- Blue vinyl (waterproof seams)
- Poly-Guard™ Special bonded waterproof coating (Super-support foam™ only)
- Black Conductive vinyl
- Coated Pressure Displacement Foam (Water resistant and cleanable)

### 9. Options

**System Options:**

- Velcro® type fasteners may be used to attach pads securely to table top

- 
- Straps for safety and restraint

Customer signature

**ORDER NO. TP-999**

Note: Custom products are quoted upon request and are not returnable. We encourage you to call to discuss your specific need and make full use of our technical assistance before ordering.

### Two Year Table Pad Warranty

AADCO Medical Inc., warrants the Table Pad products manufactured under its label to be free of defects in material and workmanship for two years from the date of purchase. Repair or replacement of any table pad product will be done at the sole discretion of AADCO Medical Inc., upon a no-charge basis and subject only to normal shipping and handling charges. This warranty assumes normal wear and tear and does not cover any product which has been abused or misused. This warranty represents the only warranty provided by AADCO Medical Inc., for its table pad products. AADCO Medical Inc., assumes no liability other than that expressly set forth in the foregoing warranty. No other warranty, whether express or implied, is provided. In no case shall the liability of AADCO Medical Inc., exceed the purchase or replacement price of the product.