

TABLE PADS, ARMBOARDS & ACCESSORIES



AADCO Table & Armboard Pads represent the finest product of their type available in the industry today. Offered in a variety of pad filler materials and upholstery/covering materials, which are explained and compared for your convenience below. All Advanced Table Pads help reduce patient discomfort and ischemia, and have a limited warranty against defects in material and workmanship for a period of two years from the date of purchase.

Five Types of Advanced Table Pad Materials

Lo-Sorb™ . . . Low X-ray Absorption Table Pads

When your application is imaging of vessels and soft tissue using contrast media designed for low to mid-diagnostic energy ranges, table pads can make a big difference in image quality. Lo-Sorb™ Table Pads provide the highly desirable low attenuation properties that maximize imaging results. Lo-Sorb™ Table Pads provide attenuation .065 mm Al. equivalent per inch* of pad thickness. Lo-Sorb™ Table Pads come with a standard black urethane coated nylon upholstery. Lo-Sorb™ Low X-ray Absorption Table Pads are standard on some of America's Top Original Equipment Manufacturer's Angiographic and X-ray Tables.

Pressure Displacing Foam™

Our standard material of high quality 2lb. density foam exceeds original manufacturers pad quality, with characteristics suitable for good pressure displacement and long procedures. 2 inches of Pressure Displacing Foam™ has X-ray attenuation properties equivalent to 0.2mm of aluminum. X-ray radiolucent for artifact free results.

Contour Foam™

A premium material offering the excellent body support and comfort. The technology for this type of foam was originally developed for NASA and used as seating material to relieve discomfort experienced by astronauts during extended flight space missions. 2 inches of Contour Foam has X-ray attenuation properties equivalent to 2.66mm of aluminum. X-ray radiolucent for artifact free results.

Lami-Foam™

When long procedures cause bottoming out, which occurs when pressure points on the patients body begin to feel the table top, Lami-Foam is the answer. Made from a laminated combination of Contour Foam and Super-Support foam, Lami-Foam provides a superior combination of support and comfort. Two inches of Lami-Foam has X-ray attenuation properties equivalent to 2.mm of aluminum. X-ray radiolucent for artifact free results.

Super-Support Foam™

Offering extra firm support for all procedures and conditions in which firm patient support is advantageous. This is the only material available with our waterproof Polyguard™ coating. 2 inches of Super-Support Foam™ has X-ray attenuation properties equivalent to .86 mm of aluminum. Radiolucent artifact free results. Only available as a Custom Order.

Four Types of Advanced Table Pad Upholstery

Black Nylon Lightweight and radiolucent, this tough durable material stands years of hard use and exceeds the durability and strength of our best vinyl upholstery. Its special heavy duty coating is waterproof and resists staining. Seams are sewn and water resistant. Available with all foam filler materials.

Blue Vinyl Our vinyl fabric is both durable and waterproof. Seams are heat sealed and virtually waterproof. Available only with Pressure Displacement Foam™ filler material.

Black Conductive This vinyl coated nylon reinforced material is specially manufactured to be electro-conductive. Seams are sewn and water resistant. Available in all foam filler materials.

PolyGuard™ Coating This coating technique renders a table pad impervious to all fluids including betadine. It is only available on our Super-Support Foam™ filler.

Ordering Table & Armboard Pads

AADCO Medical offers a wide variety of standard size table pads. We also offer the service of making custom designed table pads and table pad systems.

Many options are available for seam types, filler material, covering material, safety straps and Velcro® type fasteners. To the right is a list of information necessary when specifying a custom made pad.

1. Table Type & Model No _____
2. Diagram of Pad
3. Material Choice

Filler Materials: (check one)

- Contour Foam™ Available in 1.5" and 2" thicknesses.
- Super-support™ foam (must be laminated in thicknesses over 1" up to 2")
- Pressure Displacement Foam™ Available in thicknesses from 1" to 8"
- Lami-Foam™ available in 1.5", 2" and 2.5" thicknesses
- Lo-Sorb Foam™ available in 1", 1.5", 2" thicknesses

4. Upholstery Choice

Covering Materials: (check one)

- Black nylon-waterproof - highly durable (water resistant seams)
- Blue vinyl (waterproof seams)
- Poly-Guard™ Special bonded waterproof coating (Super-support foam™ only)
- Black Conductive vinyl

5. Velcro® type fasteners may be used to attach pads securely to table top, size & type.