

# ***Innovative Ergonomic Solutions For Imaging ...***

## ***Why AADCO?***

***The question on the minds of those responsible for any healthcare purchase is “Why?” Why, this company & Why this product . . .***

***For 25 years AADCO has Manufactured the highest quality Imaging Accessory Products in the Industry . . .***

***We appreciate this opportunity to answer the question,***


## ***Why AADCO?***



***Robert L. Marchione***  
***President & CEO***

***Our core belief is in partnership with our customers. It's a philosophy that runs deep at AADCO. We pride ourselves on our innovation and our ability to meet our customer's needs, by providing value added solutions for basic & complex problems. It is our fondest wish to say, “Yes We Can...” to our customers & then make it happen . . .***

***Doing business with AADCO Medical is a positive experience. I am confident that when you choose AADCO, you will find it to be so . . .***



# Why AADCO ???

## Protective Wear

*With AADCO you're dealing directly with the manufacturer. This results in a better value for you & first hand attention from those who understand the products best . . .*

*AADCO is the only manufacturer worldwide to stand behind its products with a Lifetime Warranty . . .*

*AADCO Creates and Innovates with patented products like ErgoLite® & Innovative Solutions like our Cool Comfort Vest™ . . .*

*This same thoughtful innovation and attention to detail is present throughout our products . . .*

*Fast Delivery of Products, even those that are Custom Made . . .*

*Customer Service that's focused on Giving You the Best Possible Experience . . .*

*Free Disposal by Environmentally Friendly Recycling of Lead Aprons & other attenuating materials for our customers . . .*

*Manufactured in The USA . . .*



*The Patented ErgoLite® Apron with its internal support system that keeps the weight and danger of repetitive stress from your shoulders and back...*

*It's Almost This Good...*

### **A Lifetime Warranty Nobody Else Will Give You . . .**

AADCO Medical, Inc., warrants the X-ray protective wear manufactured under its label to be free of defects in material and workmanship for the life of the product. Repair or replacement will be done at the sole discretion of AADCO Medical, Inc., upon a no-charge basis and subject only to normal shipping and handling charges. This warranty assumes normal wear and tear and does not cover any product which has been abused or misused. This warranty represents the only warranty provided by AADCO Medical, Inc., for its X-ray protective wear. AADCO Medical, Inc., assumes no liability other than that expressly set forth in the foregoing warranty. No other warranty, whether express or implied, is provided. In no case shall the liability of AADCO Medical, Inc., exceed the purchase or replacement price of the product.



*Cool Comfort™ Under Your X-ray Protective Wear . . .*

*Our Cool Comfort Vest™ does it . . .*



## **FREE... DISPOSAL OF X-RAY APRONS & SHIELDS**

Since 1994, AADCO has offered its customers responsible environmentally sound disposal of all radiation attenuating materials. This includes, Aprons, Vests, Kilts, Thyroid Protectors & all X-ray Protective Shielding Garments, as well as Drapes & Shields. The program also accepts lead glass and lead acrylic. All Metals are Recycled....

This is done without charge to our customers & without any reduction of discount on new items being purchased.

**Green...  Recycle...  Green...**

Putting any X-ray attenuating metals into landfills causes more environmentally destructive mining and processing to replace them... ..... Always Re-cycle .....

**RayShield®  
Overhead  
X-ray Barriers  
& Undertable  
Shields**



## **Why AADCO ???**

### **X-ray Barriers, Shields & Suspension Systems**

*Dealing directly with the manufacturer results in a better value for you . . .*

*We stand behind our X-ray Barriers & Suspension Systems with a Lifetime Warranty . . .*

*AADCO Creates and Innovates with patented products like MRI-Safe® Suspension Systems & Innovative Solutions Like our X-Drape® Sterile Disposable Procedure X-ray Shield . . .*

*This same thoughtful innovation and attention to detail is present throughout our products . . .*

*Fast Delivery of Products, even those that are Custom Made . . .*

*Customer Service that's focused on Giving You the Best Possible Experience . . .*

*Free Disposal by Environmentally Friendly Recycling of Lead Acrylic & other attenuating materials for our customers . . .*

*Manufactured in The USA . . .*

### **A Lifetime Warranty Nobody Else Will Give You . . .**

AADCO Medical, Inc., warrants the Overhead Barrier & Equipment Suspension Systems manufactured under it's RayShield® trademark to be free of defects in material and workmanship for the life of the product. Repair or replacement of equipment parts will be done at the sole discretion of AADCO Medical, Inc, upon a no-charge basis and subject only to normal shipping and handling charges. This warranty assumes normal wear and tear and does not cover any product which has been improperly installed, abused or misused. This warranty represents the only warranty provided by AADCO Medical, Inc, for these products. AADCO Medical, Inc., assumes no liability other than that expressly set forth in the foregoing warranty. No other warranty, whether express or implied, is provided. In no case shall the liability of AADCO Medical, Inc., exceed the purchase or replacement price of the product.



**A Broad Range of Rayshield® Mobile Barriers includes Custom Sizes . . .**

### **Sterile Disposable X-ray Drapes Eliminate up to 98% of Scattered Radiation**



**D-110 X-Drape™ 12"x17" sterile disposable X-ray Shield 0.25 mm lead equivalent with cutout. Comes with 36"x36" sterile overwrap.**

Exposure to scattered radiation is directly linked to long term health problems and is clearly one of the most serious occupational hazards, that those using X-ray Fluoroscopy face each day. The closer you are to the primary beam, the closer you are to the source of scatter radiation. Alarming, an increasing number of procedures require the Clinician to be in extremely close proximity to the primary beam, causing exposure levels to exceed safe practice guidelines. Reducing your exposure to Scatter Radiation is of the utmost importance, but how? Personal shielding such as leaded glasses and eyewear certainly help, but what about additional shielding?

AADCO Medical has solved the problem of protecting the Clinician who is required to work in close proximity to the source of Scatter Radiation. Our new RayShield® X-Drape combines full .25mm lead equivalent protection in a special soft fiber sterile drape material that can be placed directly in the sterile field between the clinician and the source of scattered radiation, even when the clinician must be in very close proximity to perform the procedure.

X-Drape® can be purchased lead free, so it can be disposed of in the regular

trash or with lead at a lower price that can be re-cycled by returning to AADCO. X-Drape® can help reduce your radiation exposure in a wide range of clinical procedures and is manufactured in special designs to fit any application.



**The Only MRI-Safe®  
Overhead Suspension  
shown with MRI-Safe™  
Lamp, Injector & Dual LCD  
Mount...**  
A wide Variety of MRI-  
Safe® Accessories are  
available...



**The Only  
Compact  
80,000 Lux  
Examination  
Lamp  
available MRI-  
Safe®...**

**Overhead LCD  
Mounts Single,  
Dual, Triple &  
Quad, available  
MRI-Safe®...**

**The Only MRI-Safe®  
Wall Mounted  
Suspension shown  
with Liebel/Tyco  
Injector Mounted...**

# **A Wide Range of AADCO Manufactured Solutions...**

**C-Gel™ Pads that  
prevent pressure  
sores during long  
procedures...**

**The Only  
Image  
Intensifier  
mounted  
X-ray  
Barrier...**

**ErgoLite® The Only  
X-ray Protective  
Garment with a  
built-in Ergonomic  
Support System...**

**The Only Mobile  
Platform X-ray  
Barrier for  
Interventional  
Procedures...**

**X-Drape® our sterile  
disposable X-ray Shield  
with circular fenestration.**

**The Dual Protector™  
The Only X-ray Protective  
Eyewear with Plastic Safety  
Lenses...**

## **Innovative Ergonomic Solutions for Imaging...**

**Certificates: ISO 13485:2003, ISO 9001:2000 Quality System; UL, CE, CB Certifications, Registered FDA, SFDA, JMH**

**Trademarks: RayShield®, ErgoLite®, X-Ban®, X-Drape® and MRI-Safe®**

**Sales Offices In the USA 800-225-9014 In Asia... Hong Kong 852-2543-4286 Outside the USA 001-802-728-3400  
for more information [www.aadcomed.com](http://www.aadcomed.com)**